


- Mesh Nebuliser**
- 1 Nebulising head
 - 2 Electrodes
 - 3 Main unit
 - 4 Indicator light
 - 5 «PRESS» button
 - 6 ON/OFF button
 - 7 Batteries Compartment
 - 8 Mains Adapter Socket
 - 9 Mouthpiece
 - 10 Adult face mask
 - 11 Child face mask
 - 12 Inserting the nebulising head
 - 13 Opening the nebulising head
 - 14 Filling the nebulising head
 - 15 Connecting the mouthpiece
 - 16 Connecting the mask
 - 17 Removing the nebulising head

Dear Customer,
 This mesh nebuliser is intended for respiratory therapy for asthma, chronic bronchitis, and other inhalatory illnesses. If you have any questions or need spare parts please contact your local Microlife-Customer Service. Your dealer or pharmacy will be able to give you the address of the Microlife dealer in your country. Alternatively, visit the internet at www.microlife.com where you will find a wealth of invaluable information on our products. Retain instructions in a safe place for future reference. Stay healthy – Microlife AG!

- Read the instructions carefully before using this device.
- 1. Important Safety Instructions**
- Follow instructions for use. This document provides important product operation and safety information regarding this device. Please read this document thoroughly before using the device and keep for future reference.
 - This device may only be used for the purposes described in these instructions. The manufacturer cannot be held liable for damage caused by incorrect application.
 - Never touch the metal mesh in the nebulising head with your hands or foreign objects.
 - This device should only be used with original accessories as shown in these instructions.
 - This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the Technical Specifications section.
 - Only use the medication prescribed for you by your doctor and follow your doctor's instructions with regard to dosage, duration and frequency of the treatment plan.
 - The liquid for nebulising should be prepared based on 0.9% saline or whatever advised by your doctor. No tap or boiled water can be used. The container used for the preparation of liquid must be sterilized.
 - Never use liquids which contain ester, fatty or suspended particles (including oil/infusion of herbs); specifications may vary with the medication type of the therapy.
 - If the device is not going to be used for a prolonged period the batteries should be removed.
 - Do not use this device if you think it is damaged or notice anything unusual.
 - Never open up the device.
 - Protect it from:
 - water and moisture
 - extreme temperatures
 - impact and dropping
 - contamination and dust
 - direct sunlight
 - heat and cold
 - flammable gas, heating units and open fire
 - corrosive liquids

- To avoid any risk of electric shock when using this device please observe the following:
 - Never use the device in a humid environment such as a bathroom.
 - Never plug in the power cord or the main unit has gotten wet.
 - Never plug in or unplug the adapter cable with wet hands.
 - Check that the supply voltage at the power socket is correct, as it can cause overload to the device.
 - Unplug the device from the power supply before cleaning.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations. Keep a minimum distance of 3.3 m from such devices when using this device.
- Ensure that children do not use this device unsupervised; some parts are small enough to be swallowed. Be aware of the risk of strangulation in case this device is supplied with cables or tubes.

- 2. Preparation before Use**
- Inserting the batteries (2)
 After you have unpacked your device, first insert the batteries. The battery compartment (2) is on the bottom of the device. Insert the batteries (2 x 1.5 V, size AA), using the indicated polarity.
- Using a Mains Adapter
 Connect the device using the Microlife mains adapter (DC 3V, 1.33A).
- Only use the Microlife mains adapter available as an original accessory appropriate for your supply voltage.
 - Ensure that neither the mains adapter nor the cable are damaged.
 - Plug the adapter cable into the mains adapter socket (8) in the main unit (3).
 - Plug the adapter plug into the wall socket.

- 3. Using the Device for the First Time**
- Prior to using the device for the first time, we recommend cleaning it as described in the section 'Cleaning and Disinfecting'.
- 1 Snap the nebuliser lock (1-B) upward to open the nebulising head (3).
 - 2 Fill the medication as per your doctor's instructions into the drug storage area (14). Ensure that you do not exceed the maximum fill volume.
 - 3 Press the nebuliser lock (1-B) downwards; a «click» sound indicates that the nebuliser cover (1-A) is closed correctly.
 - 4 Connect the mouthpiece or a face mask as shown in picture 19 and 16.
 - 5 The mouthpiece gives you a better drug delivery (10) to your nostrils.
 - 6 Choose between adult (3) or child face mask (11) and make sure that it encloses the mouth and nose area completely.
 - 7 Start the nebulisation by pushing the ON/OFF button (6); the indicator light (4) will light up green.
 - 8 Breathe in and out calmly during the therapy. Sit in a relaxed position with the upper body upright. Do not lie down while inhaling. Stop inhalation if you feel unwell.
 - 9 To stop the nebulisation during the therapy, press the ON/OFF button (6); the indicator light will automatically shut down, when the liquid has been fully used.
 - 10 If the Microlife mains adapter is used, unplug the cable from the mains adapter socket (8).
 - 11 Remove the remaining medication from the nebuliser and clean the device as described in the section 'Cleaning and Disinfecting'.
 - 12 Press the «PRESS» button (5) on the rear of the main unit (3) and gently pull the nebulising head (1); out (17).
 - 13 Empty the remaining medication from the nebuliser and clean the device as described in the section 'Cleaning and Disinfecting'.

Guarantee Card
 Name of Purchaser / Nom des acheteurs / Nome del comprador / Nom de Vendeur / Alchimist adı / Nombre del comprador / Nom do comprador / Naam koper / Ονομαστικός αγοραστής / اسم المشتري

Serial Number / Numéro de série / Numero di serie / Serien-Nr. / Sen numari / Número de serie / Seriennummer / Серийный номер / آي دي / آي دي / آي دي

Date of Purchase / Date d'achat / Data de achizitie / Kaufdatum / Datum van aankoop / Fecha de compra / Datum da compra / Seriennummer / آپدیت / آئی ڈی / آئی ڈی / آئی ڈی

Specialist Dealer / Revendeur / Categoria rivenditore / Fachhändler / Uzman satıcı / Distributore specializzato / Revendedor autorizado / Spécialistes Dealer / Εξουσιοδοτημένος πωληστής / التاجر المصرح به

Pu Yuan Biotech Co., Ltd.
 4F., 4F.-1, No. 51, Ln. 33, Jhu Rd., Nanhui Dist., Taicang, Tawain, R.O.C.

EN
 Do not shake the device during operation. This may cause the device to malfunction and/or shut down.

- 4. Cleaning and Disinfecting**
- Thoroughly clean all components to remove medication residuals and possible impurities after each treatment.
- 1 Always wash your hands well before cleaning and disinfecting the accessories.
 - 2 Clean the main unit with a soft, dry cloth. Do not expose the main unit to water or heat.
- Cleaning with water**
- 1 Snap the nebuliser lock (1-B) upward to open the nebulising head (13). Fill the drug storage area with hot tap water (max. 70 °C). Close the nebulising head and shake for 30 seconds.
 - 2 Run the device for 12 minutes in order to clean the nebulising head.
- Wash all nebuliser components** (except main unit) **separately** using warm tap water (max. 70 °C) for 30 seconds. Add 5 minutes adding if necessary a small quantity of detergent following dosage and use limitations as provided by detergent manufacturer.

- Turbo cleaning function**
- We recommend using the Microlife mains adapter, as the power consumption of the turbo cleaning function is very high.
- 1 Snap the nebuliser lock (1-B) upward to open the nebulising head (13). Fill the drug storage area with hot tap water (max. 70 °C). Close the nebulising head and shake for 30 seconds.
 - 2 Plug the adapter cable into the mains adapter socket (8) in the main unit (3).
 - 3 Plug the adapter plug into the wall socket.
 - 4 Press the ON/OFF button (6) for at least 5 seconds until the indicator light (4) starts blinking green.
 - 5 Run the turbo cleaning function for 15 minutes. The device will shut down automatically.

- Disinfecting by using disinfectants**
- All nebuliser components (except main unit) can be disinfected with chemical disinfectants following dosage and use limitations as provided by disinfectant manufacturer.
- 5. Battery Replacement**
- 1 Open the battery compartment (2) on the bottom of the device.
 - 2 Replace the batteries – ensure correct polarity as indicated by the symbols in the compartment.
 - 3 Batteries and electronic devices must be disposed of in accordance with the locally applicable regulations, not with domestic waste.
- 6. Maintenance, Care, and Service**
- We recommend replacing the nebulising head, face masks and mouthpiece every year (depending on the usage).

7. Malfunctions and Actions to take

Description	Symptoms/Possible causes	Solutions
Lowering light	1. Indicator light is not lit.	3. Replace the batteries correctly.
Device not working	2. Nebulising head is blocked.	Clean the nebulising head - see section 4 'Cleaning and Disinfecting'.
No power	1. Batteries polarity are reversed.	1. Reinstall batteries correctly.
No power	2. Low battery voltage.	2. Replace the batteries.
No power	3. Adapter is not correctly connected to the main unit or the power socket.	3. Reconnect the cable correctly to the main unit and the power socket.
No power	4. Incorrect installation of nebulising head.	4. Reinstall the nebulising head correctly following the instructions (Preparation before Use).
No power	5. Incorrect installation of the mains adapter.	5. Reinstall the mains adapter correctly following the instructions (Preparation before Use).
Indicator light flashing green	1. Contamination to the nebulising head, replace if necessary.	1. Reinstall the nebulising head correctly following the instructions (Preparation before Use).
Indicator light flashing red	2. Contamination to the nebulising head, replace if necessary.	2. Clean the nebulising head, replace if necessary.
Indicator light flashing yellow	3. The liquid is not in work contact with the nebulising head.	3. Gently shake the liquid together at the nebulising head.
Indicator light flashing blue	4. Gently shake the liquid together at the nebulising head.	4. Gently shake the liquid together at the nebulising head.
Indicator light flashing purple	5. Dirty electrode pole of the main unit.	5. Clean the electrodes of the main unit with a soft, dry cloth.

- 8. Guarantee**
- This device is covered by a 2 year guarantee from the date of purchase at the place of purchase, at our discretion. Microlife will repair or replace the defective product free of charge.
- 9. Technical Specifications**
- Nebulisation rate:** 0.2 ml/min
 Particle size: = 5 µm (MMAD)
Acoustic frequency: 103 - 123 kHz ± 1k
Acoustic noise level: < 30 dB
Nebulisation capacity: min. 2 ml, max. 8 ml
Operating conditions: 10 - 40 °C / 50 - 104 °F
 30 - 85% relative humidity
Storage conditions: -20 - 70 °C / -4 - 158 °F
 30 - 85% relative humidity
Weight: 100 g (net batteries)
Dimensions: 121 x 44 x 59 mm
Automatic switch-off: Automatically shut down, when the liquid has been used up.
Voltage usage: 2 x 1.5 V alkaline batteries; size AA
Mains adapter (DC 3V, 1.33A)
 approx. 90 minutes (using new batteries)
IP Class: IP22
Reference to standards: EN 13544-1; EN 60601-1; EN 60601-1-2; ISO 10933-5; ISO CE 1639
Expected service life: Device: 2 years
 Nebulising head: 12 months
 When using the remaining medication for 15-20 minutes each time
 Technical alterations reserved.

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